

BARRIER® Clean Air Suit

Proven to maintain ultra clean air in the OR

BARRIER® Clean Air Suit can maintain ultra clean air in the operating room¹⁻⁴ (OR). It's a single-use clothing system proven to have a better protective effect than reusable systems made of Mertex, olefin, cotton or polyester^{1,4,5}.



We've heard how decreasing counts of airborne bacteria in the OR plays an important part in infection control. And we understand that people in the OR are a main source of airborne microorganisms⁶ that can cause surgical site infections. Our BARRIER Clean Air Suit is proven to contribute to hospital infection control – the below evidence summaries explain how.

Kasina P et al, 2016¹

Prospective study/clinical setting

Key findings:

- The BARRIER Clean Air Suit achieved a significantly larger reduction in colony forming units (CFU/m³) compared to the olefin suit, and also compared to the laminar air flow-assisted Mertex suit (p=0.009)

Objective:

To investigate if BARRIER Clean Air Suit alone can achieve similar air quality in the OR to reusable clothing systems assisted by mobile laminar air flow.

Methods:

- Setting: OR (hip and knee arthroplasties)
- No. of procedures: 37
- Outcome measure: CFU/m³ of air (mean, median, range)

Products tested:

- BARRIER Clean Air Suit (n=11)
- Olefin suit (n=13): reusable suit made out of olefin fabric (woven polypropylene). Test was done without laminar air flow devices
- Mertex P-3477 suit (n=13): reusable suit made out of mixed material (cotton/polyester). Test was done with using 2 TOUL mobile laminar air flow devices

Tammelin A et al, 2013⁵

Prospective studies/clinical and laboratory settings

Key findings:

- BARRIER Clean Air Suit maintained significantly lower counts of CFU/m³ compared to a reusable clothing system (Mertex) both in the OR and in an air chamber (p<0.05)
- BARRIER Clean Air Suit improved air quality (<10 CFU/m³) even when ordinary ventilation was used in the OR
- The results show that clean air suits of different materials meeting the requirements of the EN 13795 standard may still exhibit significant differences in their protection capacity (bacterial penetration)
- BARRIER Clean Air Suit has a better protective effect than a reusable system made of cotton/polyester

Objective:

To investigate if there is a difference in protective efficacy between BARRIER Clean Air Suit and a reusable clothing system.

Methods:

- Settings: OR/dispersal chamber
- No. of procedures: OR, n=10; dispersal chamber, n=5
- Outcome measure: CFU/m³ of air (mean, range)

Products tested:

- BARRIER Clean Air Suit
- Mertex P-3477 suit – reusable, mixed material (cotton/polyester/carbon fibre)

Ljungqvist B and Reinmuller B, 2012³

Laboratory study

Key findings:

- The BARRIER Clean Air Suits all resulted in a >75% reduction in CFU concentration and a >90% reduction in particles emitted from the test subjects, compared to the mixed cotton/polyester reusable material (Mertex P-3477)
- When low concentrations of airborne bacteria-carrying particles are necessary to avoid hospital infections to patients, surgical clothing systems of evaluated disposable non-woven material is preferable

Objective:

To evaluate and compare the protective efficiency of a single-use and a reusable clothing system in a laboratory setting.

Methods:

- Setting: dispersal chamber
- No. of procedures: 5 per product (4 measurements per individual)
- Outcome measures: source strength – total particulates (numbers); bacteria-carrying particles (CFU)

Products tested:

- Older version of BARRIER Clean Air Suit (shirt tucked in)
- New BARRIER Clean Air Suit, antistatic-treated (shirt outside pants)
- New BARRIER Clean Air Suit, not antistatic-treated (shirt outside pants)
- Reusable clothing system (cotton/polyester): Mertex P-3477 clean air suit. Reusable, laundered up to 50 times (shirt tucked in)

Ljungqvist B and Reinmuller B, 2013⁴

Laborative study/clinical setting

Key findings:

- BARRIER Clean Air Suits achieved a CFU concentration (CFU/m³, mean (SD)) of 11.5 ± 6 and a source strength of 1.15 ± 0.6 under the stated conditions
- The BARRIER Clean Air Suit showed higher protection efficiency when compared to that of a suit from reusable material, Mertex P-3477 (data from a previous study in the same environment)

Objective:

To evaluate the protective efficiency of a single-use clothing system in a clinical setting.

Methods:

- Setting: OR (hip and knee arthroplasties)
- No. of procedures: 10
- Outcome measures: source strength – total particulates (numbers); bacteria-carrying particles (CFU)

Products tested:

- BARRIER Clean Air Suit

Definitions

Colony forming units (CFUs): Number of particles that carry bacteria and can give rise to a colony on a culture plate.

Laminar air flow device: Instrument that creates a streamlined air flow in the operating room. This makes all particles move in the same direction.

Dispersal chamber: Test chamber with HEPA-filtered air supply and with exhaust air in which the concentration of the total number of particles and bacteria-carrying particles from the test subjects are measured in order to calculate the source strength.

Source strength: The average number of CFUs released per second from one person wearing a specified clothing system.

What is needed to limit the levels of air-borne bacteria?

- Staff clothing that reduces the emission of skin scales²
- Adequate hygienic standards⁷
- Good ventilation⁷



BARRIER® Clean Air Suit

Proven to maintain ultra clean air in the OR

References:

1. Kasina, P., Tammelin, A., Blomfeldt, A.-M., Ljungqvist, B., Reinmuller, B., Ottosson, C. Comparison of three distinct clean air suits to decrease the bacterial load in the operating room: an observational study. *Patient Safety in Surgery* 2016;10(1): DOI 10.1186/s13087-015-0091-4. 2. Tammelin, A., Hambraeus, A., Stahle, E. Routes and sources of *Staphylococcus aureus* transmitted to the surgical wound during cardiothoracic surgery: possibility of preventing wound contamination by use of special scrub suits. *Infection Control and Hospital Epidemiology* 2001;22(6):338–346. 3. Ljungqvist, B., Reinmüller, B. PEOPLE AS A CONTAMINATION SOURCE. Surgical clothing systems for operating rooms – a comparison between disposable non-woven and reusable mixed material. Technical Report from Chalmers University of Technology, Gothenburg, Sweden. 2012. 4. Ljungqvist, B., Reinmüller, B. PEOPLE AS A CONTAMINATION SOURCE. Performance of single-use surgical clothing systems for operating rooms, Technical Report from Chalmers University of Technology, Gothenburg, Sweden. 2013. 5. Tammelin, A., Ljungqvist, B., Reinmüller, B. Single-use surgical clothing system for reduction of airborne bacteria in the operating room. *Journal of Hospital Infection* 2013;84(3):245–247. 6. Hambraeus, A. Aerobiology in the operating room – a review. *Journal of Hospital Infection* 1988;11(Supplement A):68–76. 7. Hirsch, T., Hubert, H., Fischer, S., Lahmer, A., Lehnhardt, M., Steinau, H.U., Steinstraesser, L., Seipp, H.M. Bacterial burden in the operating room: impact of airflow systems. *American Journal of Infection Control* 2012;40(7):E228–E232.

These synopses of published articles have been compiled by Mölnlycke as a service to healthcare professionals. It does not contain the complete text and Mölnlycke makes no representation as to its completeness in addressing all issues in the articles.

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00. The Mölnlycke and BARRIER trademarks, names and logos are registered globally to one or more of the Mölnlycke Health Care Group of Companies. © 2017 Mölnlycke Health Care AB. All rights reserved. HQSU0205

