

PRODUCT DATA SHEET
MESTOPORE®

Product name	Mestopore®				
Product description	Mestopore is a non-sterile highly absorbent non-irritant, breathable stomacap with a fluid barrier for continent stoma. Do not dispose in WC.				
Intended use	Self-adhesive absorbent dressing for continent stomas				
Product Codes	Art.No	Product Name	Size	Wound Pad	Pack count CON/RET/TRP
	665000-03	Mestopore 9x10 cm	9x10 cm	4.5x5.5 cm	20/20/160
Manufacturer	Mölnlycke Health Care AB Box 13080 SE-40252 Göteborg Sweden				
Country of origin	Finland				
Product classification according to MDD	Class IIa				
CE-marked	Yes				
Products material	Backing	Polyurethane film laminated with polyester nonwoven			
	Wound Pad	Polypropylene nonwoven, superabsorbent fibres/pulp/polyethylene and polyester fibres			
	Wound contact layer	Polyamide net			
	Adhesive	Water based acrylic adhesive			
	Release film	Silicone coated paper			
Shelf life	N/A				
Physical /Chemical characteristics	Mestopore is released in accordance to the following standards				
	Quality Property	Reference Method	Value / Result		
	Barrier to microorganisms	ASTM F 1671	Breathable film approved		
	Adhesion to glass	T-1082	1,7 N		
	Free swell absorptive capacity	EN 13726-1 3.2	64,0 g/100cm ²		
	Fluid handling capacity	EN 13726-1 3.3	Moisture vapour transmission rate	2,0 g/10cm ² /24h	
			Fluid absorbed	4,0 g/10cm ² /24h	
Fluid handling capacity			6,0 g/10cm ² /24h		

This document has been printed by the PRIME system.
 The validity of this document cannot be guaranteed.

© Mölnlycke Health Care AB. This document is the property of Mölnlycke Health Care and must not be reproduced, disclosed to any third party or used in any unauthorised manner without written consent.

Title: PDS Mestopore

Page 2(2)

Packaging	The purpose of the shelf container and transport box is to protect material from dust, moisture and damage during transport, handling and storage.
Storage Condition	N/A

GS1 Codes	Art No	CON	RET	TRP
	665000-03	7332551866563	7332551866563	7332551866556